


TD-3140

UPPER ARM BLOOD PRESSURE MONITOR

Mb^H
Mobile Health



Irregular Heart Beats Detection



HI

Fast Heart Beats
Heart rate: > 110 bpm.



LO

Slow Heart Beats
Heart rate: < 50 bpm.



IHB

Irregular Heart Beats Heart rate: premature, delayed contraction, skipped, pause between heartbeats, extra systole.

SPECIFICATIONS

Accuracy of Pressure	±3 mmHg or ±2% of reading
Accuracy of Heart Rate	±4% of reading
Measurement Range	Systolic: 60 ~ 255 mmHg / Diastolic: 30 ~ 195 mmHg Pulse rate: 40~199 / minute
Memory Capacity	200 sets
Communication	Bluetooth 4.0
Power Source	4 x AAA
Dimension	145 (L) x 64.5 (W) x 29.6 (H) mm
Weight	240 g (without battery)
Operating Condition	5 °C ~ 40 °C; Humidity 15% ~ 93%
Battery Life	At least 200 measurements
Power Saving	180 seconds
Backlight	LCD backlight
Net Weight per Carton	12.63 kg
Cuff Size	Wide Range: 24 ~ 43 cm, M Size: 24 ~ 35 cm

SHIPPING INFO

Units per Carton	30
Gross Weight per Carton	15.6 kg
Carton Size	42 (L) x 39 (W) x 41 (H) cm



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.

www.taidoc.com