


4G REMOTE MONITORING DEVICES




Blood Pressure Monitor



Blood Glucose Meter



Weight Scale



4G MODULE

Part. No.	EXS62-W
Description	LTE 4G Module for Global Deployment
Frequency Bands	LTE Cat NB1/2: 600 (Bd71), 700 (Bd12, Bd13, Bd28, Bd85), 800 (Bd18, Bd19, Bd20, Bd26), 850 (Bd5), 900 (Bd8), AWS-3 (Bd66), AWS-1 (Bd4), 1800 (Bd3), 1900 (Bd2, Bd25), 2100 (Bd1) LTE Cat M1: 700 (Bd12, Bd13, Bd28, Bd85), 800 (Bd18, Bd19, Bd20, Bd26, Bd27), 850 (Bd5), 900 (Bd8), AWS-3 (Bd66), AWS-1 (Bd4), 1800 (Bd3), 1900 (Bd2, Bd25), 2100 (Bd1)



SPECIFICATIONS

Part no.	TD-4257B
Enzyme Type	GDH-FAD
Sample Size	0.5 µL
Reaction Time	5 seconds
Measurement Range	20 ~ 650 mg/dL (1.1 ~ 36.1 mmol/L)
Hematocrit Range	20% ~ 65%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	±15 mg/dL if < 100 mg/dL; ±15% if ≥ 100 mg/dL
Communication	4G
Ketone Warning	Yes
Power Source	Rechargeable Li-ion Battery
Memory Capacity	500 sets
Day Average	7, 14, 21, 28, 60, 90 days for blood glucose
Daily Alarm	4 daily alarms
Dimension	102 (L) x 64 (W) x 29.5 (H) mm
Weight	81g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), 10% ~ 85% R.H.
Storage Condition	Meter: -20°C (-4°F) ~ 60°C (140°F); Strip: 2°C (35.6°F) ~ 30°C (86°F)
Package	Vial pack / Single foil pack



SPECIFICATIONS

Part no.	TD-3261
Power Source	Rechargeable Li-Ion Battery and External DC 6V
Memory Capacity	Maximum 864 memory records
Communication	4G (TD-3261G)
Power Saving	180 seconds
Size of Monitor w/o Cuff	147mm (L) x 105mm (W) x 80mm (H), 500g (include battery)
Operating Condition	10°C ~ 40°C
Accuracy of Pressure	±4% of reading
Accuracy of Heart Rate	Systolic : 60 ~ 250 mmHg / Diastolic : 30 ~ 180 mmHg
Measurement Range	Pulse rate : 40 ~ 199 / minute



SPECIFICATIONS

Part no.	TD-2555
BMI (Body Mass Index)	Yes
Max. capacity	250 kg
Accuracy	4.0 kg ~ 60.0 kg ± 0.3 kg, 60.1 kg ~ 250.0 kg ± 0.5%
Graduation	0.1 kg
Age Range	10 to 80
Height Range	100 to 220 cm
User	5
Communication	4G
Memory	135 sets
Backlight	LCD backlight
Power Source	4 x AAA
Dimension	340(L) x 400 (W) x 26 (H) mm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.