



POCT BARCODE GLUCOSE METER



Proactive QC reminder



5" full-color touch screen



Graphical interface



Support secured HTTPS



2D Barcode scanner

SPECIFICATIONS

Meter

Display	5" full-color
Power	3.7V 1100mAh Li-ion rechargeable battery
Communication	Wi-Fi
Memory	Measurement: 1000 QC: 1000 (3 levels/record)
Sync Time	less than 1s / data
Security	HTTPS
Medical Command	Supported

Power Saving	after 10 minutes if left idle
Operating Conditions	8°C~45°C (46.4°F~113°F) R.H. 10~90%
Storage Conditions	-20°C~60°C (-4°F~140°F) R.H. 10~90%
Drop Test	>5 mins continuous test (protective case included)
Dimension	152.4 (L) x 70 (W) x 24.6 (H) mm
Weight	210.6g (including battery)

Docking Station

Part No.	TD-9019
Wi-Fi Spec.	802.11/a/b/g/n
Ethernet Port	10BASE-T/100BASE-TX IEEE 802.3
Power	Output: 5V DC 2A 10W Input: 100~240V AC, 50~60Hz
Operating condition	10°C to 40°C (50°F to 104°F), 15% to 93% R.H.
Dimension	110.5 (L) x 95.2 (W) x 96.5 (H) mm
Weight	150 g

Strip

Enzyme Type	GDH-FAD
Sample Size	0.5 µL
Reaction Time	5 sec
Measurement Range	10~800mg/dL (0.56~44.4 mmol/L)
Hematocrit Range	0~70%
Precision	SD < 5mg/dL (0.278mmol/L) at < 100 mg/dL (5.55 mmol/L) CV < 5% at ≥ 100 mg/dL (5.55 mmol/L)
Accuracy	±15 mg/dL if < 100 mg/dL; ± 15% if ≥ 100 mg/dL



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.