


TD-5052

AKTIVSAFE LANCET




FEATURES


- Easy to use
- Five different size
- Single-use only
- Lancet automatic retraction
- FDA 510(k) cleared Class II medical device



Twist off the protective cap and remove the cap



Press the safety lancet body firmly against the puncture site to activate the safety lancet



Dispose the used safety lancet and collect blood sample

Gauge	21G	23G	26G	28G	30G
Lancet color	 Green	 Blue	 Yellow	 Orange	 Purple
Penetration Depth	2.2 mm	2.2 mm	1.8 mm	1.8 mm	1.6 mm
Blood Volume	60-125 μ L	30-75 μ L	10-40 μ L	5-30 μ L	3-20 μ L
Minimum Blood Drop					
Packaging	100 pcs/box	100 pcs/box	100 pcs/box	100 pcs/box	100 pcs/box

SHIPPING INFO

Units per Carton	48
Gross Weight per Carton	14.5 kg
Carton Size	44.8 (L) x 37 (W) x 34.8 (H) cm



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.