


TD-5086 AKTIVSAFE LANCET



Easy to use



Single-use only




Five different size




Lancet automatic retraction




FDA 510(k) cleared Class II medical device



Twist off the protective cap and remove the cap



Press the safety lancet body firmly against the puncture site to activate the safety lancet



Dispose the used safety lancet and collect blood sample

GAUGE

21G

23G

26G

28G

30G



Lancet color	Green	Pink	Orange	Blue	Yellow
Penetration Depth	2.2 mm	2.2 mm	1.8 mm	1.8 mm	1.5 mm
Blood Volume	60-125 μ L	30-75 μ L	10-40 μ L	5-30 μ L	3-20 μ L
Minimum Blood Drop					
Packaging	100 pcs/box	100 pcs/box	100 pcs/box	100 pcs/box	100 pcs/box

SHIPPING INFO

Units per Carton	100
Gross Weight per Carton	10.9 kg
Carton Size	90 (L) x 63 (W) x 68 (H) cm



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.